ACCU-CHEK®

Aviva

**BLOOD GLUCOSE METER /
MEDIDOR DE GLUCEMIA**

**Owner's Booklet
Manual del usuario**

Dear ACCU-CHEK[®] System Owner,

Thank you for choosing the ACCU-CHEK Aviva system!

Congratulations on your decision to take control of your diabetes. We've designed your new ACCU-CHEK Aviva system with comfort, convenience, and control in mind – to help make living with diabetes a little easier.

This booklet will help you get the most from your ACCU-CHEK Aviva system.

To start testing quickly, you can also refer to the First Time Guide.

If you have questions, we are here to help. Just call the ACCU-CHEK Customer Care Service Center toll-free at 1-800-858-8072. We offer assistance 24 hours a day, 365 days a year in many languages. You can also visit www.accu-chek.com for diabetes management tools and product demonstrations.

Thanks again for choosing the ACCU-CHEK Aviva system.

The ACCU-CHEK[®] Aviva System

The ACCU-CHEK Aviva System is designed to quantitatively measure the concentration of blood glucose (sugar) using ACCU-CHEK Aviva Test Strips. The system is indicated for use by the professional and the lay user. Professionals may use the test strips to test capillary, venous, arterial, and neonate whole blood samples; lay use is limited to capillary whole blood testing. Capillary whole blood for testing of blood glucose (sugar) can be obtained from fingertip and approved alternative sites (e.g. forearm); refer to the Alternative Site Testing (AST) section of this booklet for approved alternative sites and associated limitations.

For use only with the ACCU-CHEK Aviva Test Strips and Control Solutions.

The system includes:

- **ACCU-CHEK Aviva meter with battery**
- **ACCU-CHEK Aviva test strips and code key***
- **ACCU-CHEK Aviva control solution***
- **ACCU-CHEK Multiclix lancing device***
(with a black cap for fingertip testing and a clear cap for non-fingertip testing)
- **ACCU-CHEK Multiclix lancet drum***

*may be sold separately

Need Help?

Just call the ACCU-CHEK Customer Care Service Center toll-free at 1-800-858-8072. We offer assistance 24 hours a day, 365 days a year in many languages. You can also visit www.accu-chek.com for diabetes management tools and product demonstrations.

Please complete your Warranty Card and mail it, so you receive the best customer service possible and product update news.

Before You Start Testing

About the meter and test strips

- Carefully read and follow the instructions in the Owner's Booklet and Package Inserts for the test strips and control solutions. It is very important to follow the instructions in order to prevent a wrong result or improper treatment.
- Set the time and date on your meter before you begin testing.
- Inspect the container of test strips before using them for the first time. If you see any damage to the container cap or if anything prevents the cap from closing properly, do not use the test strips. Contact the ACCU-CHEK Customer Care Service Center. Damaged test strips can cause inaccurate results, which could lead to improper treatment.

About your new meter

- Your new meter is designed and can be used for testing whole blood samples. It should not be used to diagnose diabetes.
- Although you always apply whole blood to the test strip, your system has been calibrated to deliver plasma-like values for easier comparison to lab results.
- The meter, test strips, and control solutions are only for use outside the body (in vitro). Do not eat the test strips. Do not swallow or inject control solutions or use control solutions for any purpose other than testing the ACCU-CHEK Aviva system.

About Testing Yourself or Others

- **Certain substances can interfere with the ACCU-CHEK Aviva blood glucose system and cause falsely high results. For example, peritoneal dialysis solutions containing icodextrin (such as Extraneal) or certain immunoglobulin therapies that contain maltose (such as Octagam 5 %), cause inaccurate results. For more information, refer to the Test Strip Package Insert, or ask your doctor.**
 - **If you are very dehydrated or urinating frequently, you may get an inaccurate test result. If you think you are dehydrated, call your doctor right away.**
-
- Some people with diabetes do not experience symptoms of low blood sugar (hypoglycemia). Others, such as children or people who are unconscious or have certain disabilities, may not be able to communicate their symptoms to caregivers. For these reasons, do not change any treatment without first talking to a doctor.
 - Run a control test when you open a new box of test strips or if you think that your test result is incorrect. Running a control test lets you know that the meter and test strips are working properly.
 - Refer to your Test Strip and Control Solution Package Inserts for additional health-related information.

- It is always a good idea to have a back-up testing method available. Failure to test could cause a delay in treatment decisions and lead to a serious medical condition. Examples of back-up testing methods include a back-up meter or testing by a lab. Ask your doctor or pharmacist about other possible back-up methods.

DO NOT CHANGE YOUR TREATMENT BASED ON A SINGLE RESULT THAT DOES NOT MATCH HOW YOU FEEL OR IF YOU BELIEVE THAT YOUR TEST RESULT COULD BE INCORRECT.

If your blood sugar result doesn't match how you feel and you have followed the instructions in this Owner's Booklet, follow your doctor's instructions, or call your doctor.

Special Information for Healthcare Providers and Caregivers

- Do not use this device to measure blood sugar in people who are experiencing cardiovascular collapse (severe shock) or decreased peripheral blood flow.
- Consult your physician to determine if it is appropriate for your child to be taught how to use the meter system or any other medical products.

Contents

Chapter 1: Understanding Your New System	9
The ACCU-CHEK Aviva Meter	9
Coding the Meter	11
Adjusting the Time and Date – First Time Use.....	13
Using the ACCU-CHEK Aviva Test Strips	15
Chapter 2: Control Testing	17
Why Run Control Tests	17
About the Control Solutions	18
Running a Control Test	19
Understanding Control Test Results	23
Chapter 3: Testing Your Blood Sugar	27
Using the ACCU-CHEK Multiclix Lancing Device.....	27
Inserting a Lancet Drum	28
Preparing the Lancing Device for Fingertip Testing.....	30
Running a Blood Sugar Test with Blood from Your Fingertip.....	32
Changing the Lancet Drum	37
Running a Blood Sugar Test with Blood from Your Palm, Forearm, Upper Arm, Thigh, or Calf	39
Flagging Test Results	48
Understanding Test Results	49
Unusual Test Results.....	49
Symptoms of High or Low Blood Sugar	50
Comparing Your Meter Result to a Lab Result	51

Chapter 4: Meter Memory, Setup, and Downloading.....	53
Memory	53
Viewing Test Results.....	55
Meter Setup.....	56
Setting the Time and Date	57
Setting the Beeper On/OFF	59
Setting Test Reminders.....	62
Setting the Hypoglycemic (Hypo) Indicator	66
Downloading Your Results to a Computer or PDA.....	69
Chapter 5: Maintenance and Troubleshooting	71
Changing the Battery	71
Cleaning the Meter	73
Maintenance and Testing.....	74
Cleaning the Lancing Device.....	74
Screen Messages and Troubleshooting	75
Chapter 6: Technical Information	83
Product Limitations.....	83
Specifications	83
Product Safety Information	85
Explanation of Symbols.....	85
Travel Documentation.....	86
Disposing of the Meter, Test Strips, Lancets, and Batteries.....	86
Warranty.....	88
Limited License	91
Additional Supplies	94
Information for Healthcare Professionals	95
Index.....	97

Chapter 1: Understanding Your New System*

The ACCU-CHEK Aviva Meter

Infrared (IR) Window – Transfers data from the meter to a computer or PDA.

Display – Shows results, messages, and results stored in memory.

Right and Left Arrow Buttons – Press to view memory, adjust settings, and scroll through results.

Test Strip Slot – Insert test strip here.

On/Off/Set Button – Turns the meter on or off and sets options.

Battery Door – Flip open the battery door by pushing the tab in the direction of the arrow.

Code Key Slot – Insert code key into this opening.

Code Key
(for example)

Battery – Insert with (+) symbol facing up.

Test Strip**Golden End –**

Insert this end of the test strip into the meter.

**Yellow Window –**

Touch blood drop or control solution here.

Test Strip Container**Control Solution Bottle****Lancet Drum****Clear Cap**

(for obtaining blood from sites other than your fingertip)

Top View

Infrared (IR) Window

On/Off/Set Button

Lancing Device

Comfort Dial with Depth Selection

Release Button

Lancet Counter

Cap

Depth Indicator

Plunger

*Some items may be sold separately.

Coding the Meter


A new code key is inside each new box of test strips.

1. Make sure the meter is off.

2. Turn over the meter.

3. Remove the old code key (if there is one in the meter) and discard it.

**Change the code key every time you open a new box of test strips!
For accurate results, it is important to be sure the code key matches
the code number on the test strip container.**

Notes:


- Don't force the code key into the meter – it is made to go into the meter only one way.
- If you see “code” and “- - -” on the display, turn off the meter and reinsert a code key into the meter.


4. Turn over the code key so the code number faces away from you. Push it into the code key slot until it stops.
5. Leave the code key in the meter until you open a new box of test strips. Remember to change the code key each time you open a new box of test strips.

Adjusting the Time and Date – First Time Use

Your new meter comes with a preset time and date. You may need to change the time to your time zone. Having the right time and date in your meter is important if you use the meter memory or if you want to download your results to a computer. It also helps your healthcare team interpret your results.

decrease

increase



1. Press and release  to turn on the meter. The time and date appear on the display. “Set-up” and the hour flash. If this is the correct time and date, press and **hold**  to exit the set-up mode.

2. Press and release  or  to decrease or increase the hour. Press and **hold**  or  to scroll faster.

3. Press and release  to set the hour. The minutes flash.

Note: When you install a new battery, the meter automatically prompts you to check the time and date when you turn it on. If the time and date are correct, press and **hold**  to exit the set-up mode.

decrease increase

4. Repeat steps 2 and 3 to set the minutes, am/pm, month, day, and year. The flashing field is the one you are changing.

5. Once you set the year, press and **hold**  until the flashing test strip symbol appears.

6. See Chapter 4, “Meter Setup” to set other options.

Using the ACCU-CHEK Aviva Test Strips

- Use only ACCU-CHEK Aviva test strips. Using other test strips with this meter can produce inaccurate results.
- Replace the code key and run a control test every time you open a new box of test strips. See Chapter 2, “Control Testing.”
- Store the unused test strips in their original container with the cap tightly closed. Do not remove test strips from the test strip container and put them into another container, such as a plastic bag, pocket, purse, wallet, etc.
- Close the container tightly right after you take out a test strip. This helps keep the test strips dry.
- Be sure to check the “Use By” date on the test strip container. Do not use the test strips after that date. See Chapter 5, “Screen Messages and Troubleshooting” for the error code for expired strip. If the “Use By” date is missing or cannot be read, do not use the test strips and call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
- Use the test strip within three minutes after you take it out of the container.
- Refer to the test strip package insert for test strip storage and system operating conditions.
- Do not apply blood or control solution to the test strip before you insert it into the meter. If the meter gives you a result before applying blood or control solution, do not act on that result.
- Do not reuse test strips. Once you have put control solution or blood on a test strip, discard it. If you need to retest, use a new test strip.

- **DO NOT** expose strips to heat, moisture or humidity. Temperatures outside the required range, as well as moisture and humidity (bathroom, kitchen, laundry room, car, or garage) can damage your test strips and lead to inaccurate results.
- **DO NOT** bend, cut, or alter the test strips.
- **DO NOT** get dirt, food, or other material on the test strip.

Not following these precautions can lead to inaccurate results.

Chapter 2: Control Testing

Why Run Control Tests

Running a control test lets you know that the meter and test strips are working properly. You should run a control test when:

- You open a new box of test strips
- You left the test strip container open or you think your test strips have been damaged
- Your test strips were stored in extreme temperatures and/or humidity
- You want to check the meter and test strips
- You dropped the meter
- Your test result does not agree with how you feel
- You want to check if you are testing correctly

About the Control Solutions

- Only use ACCU-CHEK Aviva control solutions.
- The meter is designed to recognize the difference between the ACCU-CHEK Aviva control solution and blood.
- Write the date you open the bottle on the bottle label. The control solution must be discarded three months from the date the bottle was opened (discard date) or on the “Use By” date on the bottle label, whichever comes first.
- Do not use control solution that is past the “Use By” or discard date.
- The control results are not displayed in memory.
- The control solution can stain clothing. If you spill it, wash your clothes with soap and water.
- Close the bottle tightly after use.
- Store the bottle at 36 °F to 90 °F. Do not freeze.

Running a Control Test

You need the meter, a test strip, and control solution Level 1 or Level 2. The control level is printed on the bottle label.

Note: Most people just test the Level 1 control. If you wish, you can also test a Level 2 control. A set of Level 1 and Level 2 control solutions is available for purchase. To order control solutions, talk to your pharmacist or visit www.accu-chek.com to order online. Your meter is designed to recognize the difference between the ACCU-CHEK Aviva control solution and blood. The control results are not displayed in memory.

1. To make sure the display is working properly, turn off the meter, then press and **hold**  to see the complete display. All segments should be clear and look exactly like the picture. If any segments are missing from the display, there may be a problem with the meter. Call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

2. Insert the test strip into the meter in the direction of the arrows. The meter turns on.


3. Make sure the code number on the display matches the code number on the test strip container. If you miss seeing the code number, remove the test strip and reinsert it into the meter.



4. Select the control solution you want to test. You will enter the level later in the test.


5. Place the meter on a flat surface, like a table.
6. Remove the control bottle cap and wipe the tip of the bottle with a tissue.

7. Squeeze the bottle until a tiny drop forms at the tip of the bottle. Touch the drop to the **tip** of the yellow window of the test strip. **Do not put control solution on top of the test strip.** When you see ⌚ flash, you have enough control solution in the test strip. Wipe the tip of the bottle with a tissue then cap the bottle tightly.

8. A result appears on the display, along with the control bottle symbol and a flashing “L.” **Do not remove the test strip yet.** Press and release  once to mark the result as a Level 1. If you tested the Level 2 control, press and release  a second time.

9. Press and release  to set the control level in the meter.


10. The control result and “OK” alternate on the display if the result is in range. The range is printed on the test strip container label. The control result and “Err” alternate on the display if the result is out of range. Remove and discard the used test strip.

Understanding Control Test Results

The label on your test strip container shows the acceptable ranges for both the Level 1 and the Level 2 control solutions. The result you get should be inside this range. Make sure you compare the result to the correct level of control.

When the control result is inside the range on the test strip container, the test strips and meter are working properly.

23

Example

	Range (mg/dL)
Level 1	25–55
Level 2	255–345

(This is an example. Refer to the ranges on your test strip container.)

The control ranges apply only to control results. They only indicate that your test strips and meter are working properly. Do not use control results to interpret blood sugar results.

If the control result is not inside the acceptable range (printed on the test strip container), do not use the meter until you solve the problem. Here are some things you can check:

Troubleshooting Checks**Action**

1. Are the test strips or control solution past the “Use By” date or discard date?	If either is past the “Use By” date or discard date, throw it away.
2. Did you wipe the tip of the control solution bottle with a tissue prior to use?	Wipe the tip of the bottle with a tissue. Repeat the control test with a new test strip and a fresh drop of control solution.
3. Were the test strip container and control solution bottle always closed tightly?	If you think either may have been uncapped for some time, replace the test strips or control solution.
4. Was the test strip out of the container for more than three minutes?	Repeat the control test with a new test strip.
5. Were your test strips and control solution stored in a cool, dry place?	Repeat the control test with properly stored test strips or control solution.

Troubleshooting Checks	Action
6. Did you follow all of the testing instructions?	Read Chapter 2, “Control Testing” and test again. If you still have problems, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
7. Did you choose the correct control solution level, either 1 or 2, when you ran the test?	If you chose the wrong control solution level, you can still compare the control result to the range printed on the test strip container.
8. Does the code number on the meter display match the code number on the test strip container?	If they don't match, insert the correct code key in the meter and test again.
9. If you are still unsure of the problem ...	Repeat the control test with a new test strip. If you still have problems, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

Chapter 3: Testing Your Blood Sugar

Using the ACCU-CHEK Multiclix Lancing Device

The ACCU-CHEK Multiclix lancing device is intended for patient self-monitoring by a single person only. It must not be used to collect blood from more than one person as this poses an infection risk.

Notes:

- Testing with blood from sites other than your fingertip requires the clear cap. This is explained later in the chapter.
- The best depth setting is the lowest number that lets you get enough blood for a test. Try different settings to find the one that is right for you. For virtually pain-free testing, always use a new lancet for every test. With the ACCU-CHEK Multiclix lancing device, you simply change the lancet with a twist.
- To reduce the risk of infection, never share your lancing device with anyone.
- Always use a new lancet with each test.

The black cap is for fingertip testing only.

The clear cap is for testing other approved sites on the body, not the fingertip.

A new drum has a blue ring that does not move.

A used drum has two visible red stripes and the blue ring can rotate on the white drum. You cannot insert a used drum into the lancing device.

Inserting a Lancet Drum

You must first load the drum into the lancing device to get it ready for use.

28

1. Remove the cap by pulling it straight off. It is not a twist cap. The easiest way to do this is to place your thumb to the side of the notch as seen in the illustration.
2. Before you insert a drum, be sure that the release button is **not** yellow, and that no lancets are sticking out of the drum. Insert a new drum, blue ring first, until it stops. It clicks into place.

3. Slide the cap on until it stops by aligning the notch on the cap with the notch on the device.

4. You are now ready to use the first lancet. The lancet counter on the plunger shows a number 6, meaning you have six new lancets remaining.

Preparing the Lancing Device for Fingertip Testing

Use a new lancet for each fingerstick to avoid infection and get a virtually pain-free test.


30

1. Make sure the black cap is on the lancing device. Do not use the clear cap to perform a fingertip test.
2. Adjust the depth by turning the Comfort Dial. The depth indicator shows the current depth selection. The higher the number, the deeper the penetration. If you have soft skin, we suggest starting at a depth of 2. If your skin is calloused or thick, try a higher setting.

3. Press the plunger in as far as it will go, like a pen.

4. The release button turns yellow. Do not press the release button while pressing the plunger. Set the lancing device aside until later in the test.


Running a Blood Sugar Test with Blood from Your Fingertip

Before you run your first blood sugar test, set up the meter correctly and run a control test. You need the meter, a test strip, and a lancing device with a drum loaded to run a blood test.


1. To make sure the display is working properly, turn off the meter, then press and **hold**  to see the complete display. All segments should be clear and look exactly like the picture. If any segments are missing from the display, there may be a problem
2. Wash your hands with warm soapy water and dry thoroughly. Dirty or wet hands could affect results.
3. Prepare the lancing device.
4. Insert the test strip into the meter in the direction of the arrows. The meter turns on.
5. Make sure the code number on the display matches the code number on the test strip container. If you miss seeing the code number, remove the test strip and reinsert it into the meter.

Note: After completing a blood sugar or control test, the meter turns itself off five seconds after the test strip is removed.

6. When the blood drop symbol flashes, obtain a drop of blood from your fingertip.

7. Hold the lancing device firmly against the side of your fingertip. Remember that the opening where the lancet comes out is not in the center of the cap. Press the yellow release button. Make sure that each lancet is fully retracted after use.

8. Gently squeeze your finger to assist the flow of blood. This helps you get a blood drop. Touch the blood drop to the **tip** of the yellow window of the test strip. **Do not put blood on top of the test strip.**

9. When you see  flash, you have enough blood in the test strip. If you applied blood but do not see the flashing , you may reapply more blood within five seconds.

10. The result appears on the display. If you want to flag your test result for a special event, leave the test strip in the meter. See Chapter 3, “Flagging Test Results.” Otherwise, throw the used test strip away.

If mg/dL does not appear with the test result call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072. Use of the wrong unit of measure may cause misinterpretation of your actual blood sugar level, and may lead to incorrect treatment.

Note: Do not twist the plunger when the release button is yellow or when the lancet counter shows a number 1. This may damage the device.

36

11. For virtually pain-free testing, we recommend that you advance to a new, sterile lancet after each test. To load the next lancet, twist the plunger one-quarter turn forward (until it stops), then twist it back all the way.

12. The lancet counter decreases one number. For safety reasons, once you advance to a new lancet, you cannot go back to a used lancet.

Changing the Lancet Drum

When you have used the sixth and last lancet, you should change the drum.

1. Take the black cap off the lancing device by pulling it straight off. It is not a twist cap. The easiest way to do this is to place your thumb to the side of the notch as seen in the illustration.
2. Hold the drum between your thumb and index finger and pull it straight out of the lancing device. Two red stripes are visible on the white part, indicating that the drum has been used.

Notes:

- A used drum cannot be reinserted into the lancing device.
- If the ACCU-CHEK Multiclix lancing device or drum is dropped on a hard surface, the lancets may stick out slightly from the drum and could cause injury. If you drop a drum, always pick it up at the red-striped end and dispose in a sharps container.

38


3. Throw the old drum away. Always dispose in accordance with local regulations.

4. Before you insert a drum, be sure that the release button is **not** yellow. Insert a new drum, blue ring first, until it stops. It clicks into place.

Running a Blood Sugar Test with Blood from Your Palm, Forearm, Upper Arm, Thigh, or Calf (non-fingertip testing)

You have the option of testing other sites on your body besides the fingertip. Blood obtained from a fingertip can be used at any time to measure blood sugar. The two palm testing sites are the fleshy areas under the thumb (thenar) and under the little finger (hypothear). If blood from an alternative site such as the palm, forearm, upper arm, thigh, or calf is used, there are certain times when testing is not appropriate (see next section). This is because your blood sugar level changes faster in your fingertip than in the alternative sites. These differences could cause you to make the wrong therapeutic decision, producing adverse health effects. Please read the next section before you try testing from other sites.

IMPORTANT

Talk to your doctor before you begin using other test sites.

Do not change your treatment because of just one result.

NEVER ignore symptoms of high or low blood sugar.

If your blood sugar result does not match how you feel, do a fingertip test to confirm your result. If the fingertip result still does not match how you feel, call your doctor.

Tests from sites other than your fingertip may be done:

- Immediately before a meal
- Fasting

Fingertip test only:

- Two hours or less after eating
- After exercising
- If you are sick
- If you think your blood sugar is low
- If you often don't notice when your blood sugar is low
- During peak action time of short-acting insulin or rapid-acting insulin analogues
- Up to two hours after injecting a short-acting insulin or a rapid-acting insulin analogue


You need the meter, a test strip, a lancing device with a drum loaded, and the clear cap from the kit.

1. To make sure the display is working properly, turn off the meter, then press and **hold**  to see the complete display. All segments should be clear and look exactly like the picture. If any segments are missing from the display, there may be a problem with the meter. Call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
2. Wash your hands and the test site with warm, soapy water and dry thoroughly.
3. Remove the black cap from the lancing device by pulling it straight off. It is not a twist cap. The easiest way to do this is to place your thumb to the side of the notch as seen in the illustration.
4. Put the clear cap on by lining up the notched sections on the cap and the lancing device.

Note: We suggest starting at a depth setting of 5½ for forearm, upper arm, thigh, or calf testing, and 3 for palm testing. Once you run a successful test, you may want to find the lowest depth setting that provides enough blood with the least amount of pain.

42

5. Dial the depth selection appropriate for your testing site.

6. Press the plunger in as far as it will go, like a pen. The release button turns yellow. Set the lancing device aside until later in the test.

7. Put a test strip into the meter in the direction of the arrows. The meter turns on.

8. Make sure the code number on the display matches the code number on the test strip container. If you miss seeing the code number, remove the test strip and reinsert it into the meter.
9. When the blood drop symbol flashes, obtain a blood drop from the site. Make sure that each lancet is fully retracted after use.


TIP: Rub the skin prior to lancing to increase blood flow.

44

10. Press the opening of the clear cap firmly against a fleshy area on the non-fingertip site. Press the lancing device up and down in a slow pumping motion to assist the flow of blood. Remember that the opening where the lancet comes out is not in the center of the cap.

11. Keep steady pressure on the area and press the yellow release button. Don't be afraid to apply pressure to the area with the lancing device to assist the flow of blood. Continue to press the cap against the area for a few seconds to allow blood to come to the surface. You can see the blood drop through the clear cap.

12. Touch the blood drop to the **tip** of the yellow window of the test strip. **Do not put blood on top of the test strip.**

13. When you see ⌚ flash, you have enough blood in the test strip. If you applied blood but do not see the flashing ⌚, you may reapply more blood within five seconds.

10:38^{am} 12-11


46

14. Your result appears on the display. If you want to flag your test result for a special event, leave the test strip in the meter. See Chapter 3, “Flagging Test Results.” Otherwise, discard the used test strip.
15. For virtually pain-free testing, we recommend that you advance to a new, sterile lancet after each test. To load the next lancet, twist the plunger one-quarter turn forward (until it stops), then twist it back all the way.

If mg/dL does not appear with the test result call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072. Use of the wrong unit of measure may cause misinterpretation of your actual blood sugar level, and may lead to incorrect treatment.

16. The lancet counter decreases one number. For safety reasons, once you advance to a new lancet, you cannot go back to a used lancet.

Having trouble testing from sites other than your fingertip?

If you experience difficulty obtaining a sufficient drop of blood while testing from a site other than your fingertip, please contact the ACCU-CHEK Customer Care Service Center at 1-800-858-8072 for assistance. We will guide you through the testing procedure step by step. If you continue to have difficulty, we can also offer you, free of charge, our ACCU-CHEK Softclix lancing device. Because the ACCU-CHEK Softclix lancet has a thicker diameter needle than the ACCU-CHEK Multiclix lancet, you may have better success when trying to test from sites other than your fingertip.


Flagging Test Results

You might want to flag an event such as an after-meal (post prandial) test, exercise, or a non-fingertip result. If you wish, you can “flag” your test result with an asterisk * to mark a special result. When you review your results in memory, this “flag” can help you remember what was different about the result.

48

Here is how to flag a result:

1. Run a test.
2. With the test result on the display, and the test strip **STILL IN THE METER**, press and release  or  once. The * appears in the lower right corner.
3. Remove and discard the used test strip.

Understanding Test Results

Unusual Test Results

If your blood sugar result does not match how you feel, follow these steps:

1. Run a control test. See Chapter 2, “Control Testing.”
2. Repeat a blood sugar test. See Chapter 3, “Running a Blood Sugar Test with Blood from Your Fingertip.”
3. If your blood sugar results still do not match how you feel, follow your doctor’s instructions or call your doctor immediately.

Note: Always follow your doctor’s instructions. For example, if your doctor has advised you to immediately treat a low blood sugar result (such as by eating something), then do that first.

Symptoms of High or Low Blood Sugar

Being aware of the symptoms of high or low blood sugar can help you understand your test results and decide what to do if they seem unusual. Here are some common symptoms:

High blood sugar (hyperglycemia): increased fatigue, increased appetite or increased thirst, frequent urination, blurred vision, headache, general aching, or vomiting.

Low blood sugar (hypoglycemia): sweating, trembling, blurred vision, rapid heartbeat, tingling, or numbness around mouth or fingertips.

If you are experiencing any of these symptoms, or other unusual symptoms, test your blood sugar. If your blood sugar result is displayed as LO or HI and you have symptoms of low or high blood sugar, follow your doctor's instructions or contact your doctor immediately. If your blood sugar result does not match how you feel, follow the steps in Chapter 3, "Unusual Test Results."

The meter is designed to display blood sugar results from 10 to 600 mg/dL. If a result below 10 mg/dL or over 600 mg/dL appears on the display, call the ACCU-CHEK Customer Care Service Center toll-free at 1-800-858-8072.

Comparing Your Meter Result to a Lab Result

A common question is how the blood sugar results on your meter compare to the lab results. Your blood sugar can change quickly, especially after eating, taking medication, or physical activity. If you test yourself in the morning and then go to the doctor's office for a blood sugar test, your results will probably not match even if you are fasting. This is typically not a problem with your meter; it just means that time has elapsed and your blood sugar has changed.

If you want to compare your meter result to the lab result, **you must be fasting**. Take your meter to the doctor's office, and test yourself by fingerstick within five minutes of having blood drawn from your arm by a healthcare professional. Keep in mind that the lab uses different technology than the meter, and that blood glucose meters for self-testing generally read somewhat lower than the lab result.

If you are fasting and you do a fingerstick test within five minutes of having your blood drawn, here are the general guidelines to compare your meter result to the lab result:

- If your blood sugar is below 75 mg/dL, your results generally should fall within ± 15 mg/dL of the lab result.
- If your blood sugar is equal to or over 75 mg/dL, your results generally should fall within ± 20 % of the lab result.

Chapter 4: Meter Memory, Setup, and Downloading

Memory

Storing Blood Sugar and Control Test Results

The meter automatically stores up to 500 blood sugar test results and 20 control results with the time and date of the test. You can review them at any time. Test results are stored from the newest to the oldest. It is very important to set the correct time and date in the meter. Having the correct time and date setting helps ensure appropriate interpretation of stored blood sugar results by you and your healthcare team.

Notes:


- The memory is not lost when you replace the battery. You do need to check that the time and date are correct after you replace the battery. See Chapter 4, “Setting the Time and Date.”
- Once 500 blood sugar or 20 control results are in memory, adding a new result causes the oldest one to be deleted.
- Press and **hold** ◀ or ▶ to scroll faster through the results.
- Up to 20 control results are stored in memory but cannot be reviewed on the meter. The stored control results must first be downloaded to a compatible software application. For product availability, please contact the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
- Control results are not included in the 7-, 14-, and 30-day averages.

Do not change your therapy based on an individual result in memory. Talk to your doctor before changing therapy based on results in memory.

Viewing Test Results

1. Press and release either ◀ or ▶ to access memory. The most recent result appears.
2. Press and release ◀ to view previous results in order.
3. OR, press and release ▶ to look at the 7-, 14-, or 30-day averages, in that order. You see an “n” with a number in the upper left corner of the display. This tells you how many test results are included in that average.

Meter Setup

Using the Set-Up Mode

By using the set-up mode, you can personalize the meter to suit your lifestyle. Here are the features you can customize:

Time and date – set the time and date.

Beeper – select “On” or “OFF.”

Test reminders – select “On” or “OFF.” If you choose “On,” select up to 4 times a day for a reminder to test.

Hypoglycemic (Hypo) indicator – select “On” or “OFF.” If you choose “On,” choose the blood sugar level for the indicator.

Throughout this section, the purple progress bar helps guide you through the setup process.

You are here

TIME / DATE
(hrs, min, am/pm,
month, day, year)

BEEPER
(on/off)

TEST REMINDERS
(A-1, A-2, A-3, A-4)


HYP0
(off, on, level)

Setting the Time and Date


1. Press and release  to turn on the meter. The flashing test strip symbol appears.

2. To enter the set-up mode, press and **hold**  for about four seconds. "Set-up" and the hour flash on the display.


decrease increase

- 3 Press and release  or  to decrease or increase the hour. Press and **hold**  or  to scroll faster.

4 Press and release  to set the hour. The minutes flash.

5. Repeat steps 3 and 4 to set the minutes, am/pm, month, day, and year. The flashing field is the one you are changing.

6. If you want to set up more options, press and release . If you want to exit, press and **hold**  until you see the flashing test strip symbol.

Setting the Beeper On/OFF


Your new meter has the beeper preset to “On.” You can set the beeper to “OFF,” if you prefer – this won’t affect your test results.

The beeper is helpful because it prompts you:


- To apply blood or control solution to the test strip
- When enough blood or control solution is drawn into the test strip
- When the test is complete
- When a button is pressed
- When it is time to test (if you set test reminders)
- If an error occurred while testing (even if the beeper is set to “OFF,” it still beeps for an error)

1. Press and release  to turn on the meter. The flashing test strip symbol appears.

2. To enter the set-up mode, press and **hold**  for about four seconds. “Set-up” and the hour flash on the display.

3. Press and release  repeatedly until you reach the display with the flashing beeper symbol and “On.”

4. Press and release  or  to switch between “On” and “OFF.”

5. If you want to set up more options, press and release . If you want to exit, press and **hold**  until you see the flashing test strip symbol.

TIME / DATE
(hrs, min, am/pm,
month, day, year)

BEEPER
(on/off)

TEST REMINDERS
(A-1, A-2, A-3, A-4)

HYP0
(off, on, level)

You are here

Setting Test Reminders

Test reminders are a handy way to remind you to test. You can set up to four reminders per day. The meter beeps at the time you set, again two minutes later, and two minutes after that unless you insert a test strip or press and release any button. You must have the beeper set to “On” for this option. Your meter is preset with the test reminder set to “OFF.” You must turn it on to use this feature.

If you turn A-1, A-2, A-3, and A-4 on, the meter is preset with the following times for your convenience. You can adjust the times to suit your needs.

A-1 8:00 am **A-2** 12:00 pm (noon) **A-3** 6:00 pm **A-4** 10:00 pm

Notes:


- If you tested within 30 minutes of a test reminder, the reminder does not occur.
- If the meter is on at the test reminder time, the reminder does not occur.
- When the meter beeps for a test reminder, the time, the reminder number (A-1 to A-4), and the flashing bell symbol are displayed.


1. Press and release  to turn on the meter. The flashing test strip symbol appears.


2. To enter the set-up mode, press and **hold**  for about four seconds. “Set-up” and the hour flash on the display.


3. Press and release  repeatedly until you reach the display with the bell symbol, “OFF,” and the flashing “set-up” and “A-1.” When you are setting the time for test reminders, the bell symbol remains on the display and “set-up” flashes continuously.


4. Press and release  or  to select “On” or “OFF.” Press and release  to set your choice.

5. If you select “On,” the hour flashes. “A-1” and the bell symbol remain on the display.
6. Press and release  or  to select the hour. Press and release  to set the hour.


7. The minutes flash. Press and release  or  to select 00, 15, 30, or 45. These are the only choices.

8. Press and release  to set the minutes. “Am” or “pm” flashes on the display.

9. Press and release  or  to switch between “am” and “pm.” Press and release  to set “am” or “pm.”

10. The next test reminder “A-2” and “set-up” flash on the display with “OFF” and the bell symbol. You can set additional test reminders or press and **hold**  until you see the flashing test strip symbol to exit the set-up mode.

TIME / DATE
(hrs, min, am/pm,
month, day, year)

BEEPER
(on/off)

TEST REMINDERS
(A-1, A-2, A-3, A-4)

You are here


HYPO
(off, on, level)

Setting the Hypoglycemic (Hypo) Indicator



Healthcare professionals: The hypoglycemic level may vary from person to person. It is recommended to turn **off** the hypo alarm function when using the meter in a professional facility.

66


You can set the meter to let you know when your result indicates a possible low blood sugar (hypoglycemia). You can also select what blood sugar level you want this indicator to have (60 to 80 mg/dL). **Before you set the indicator, talk to your doctor to help you decide what blood sugar level is your hypo level.** Your new meter is preset to “OFF” for the hypo indicator. If you want to turn it “On,” follow these steps.


1. Press and release  to turn on the meter. The flashing test strip symbol appears.


2. To enter the set-up mode, press and **hold**  for about four seconds. "Set-up" and the hour flash on the display.
3. Press and release  repeatedly until you reach the display with "OFF" and the flashing "set-up" and .
4. Press and release  or  to switch between "On" and "OFF." Press and release  to set your choice.

5. If you select the hypo indicator “On,” “set-up,” and  flash. The display shows 70 mg/dL.

6. Press and release  or  to select the level between 60 to 80 mg/dL. Press and release  to set it.

7. Press and **hold**  until you see the flashing test strip symbol to exit the set-up mode.

This function is no substitute for hypoglycemia training by your doctor or diabetes team.

Downloading Your Results to a Computer or PDA

We offer a variety of software to help you download your results. For information on ACCU-CHEK software, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.


You can transfer the stored results to a computer to track, identify patterns, and print.

Transferring Data Directly to a Computer or PDA Using Specialized Software and Infrared Cable

1. Install the software, if necessary, according to the instructions. To transfer the results to a computer, connect the infrared meter cable according to the instructions.
2. Run the software program and follow the instructions about how to download information. Make sure the software is ready to accept data from the meter.
3. With the meter off, press and **hold** both  and  until two arrows on the display alternately flash.
4. Locate the infrared (IR) window on the top of the meter.
5. Locate the IR window on either the infrared cable (computer) or PDA.

- Put the meter on a flat surface. Point the two IR windows toward each other. They should be 1 to 4 inches apart.
- Do not move the infrared cable (computer), PDA, or meter during the transfer.
- Follow the prompts on the software.
- The software program may shut off the meter automatically when the data transfer is complete. Should this occur, follow the prompts on your screen.

Downloading data is not recommended for healthcare facilities since there is no means to identify results by patient.

Notes:

- If the data did not transfer successfully, try again. If you still have problems, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
- To make the most of the download feature, you must set the time and date in the meter correctly.

Chapter 5: Maintenance and Troubleshooting

Changing the Battery

1. Open the battery door on the back of the meter by pushing the tab in the direction of the arrow and pulling up the door. Remove the old battery.
2. Insert the new battery with the + side up. It does not snap into place but rests on the metal contact. The door holds down the battery.
3. Put the battery door back in place and snap it closed.

Notes:

- The meter uses one 3-volt lithium battery, coin cell 2032. This type of battery can be found in many stores. Always keep a spare packaged battery on hand.
- Be sure the battery goes in + side up, or facing you.
- After you change the battery, your meter prompts you to confirm the meter's time and date settings. See Chapter 4, "Setting the Time and Date." All test results are saved in memory.

Cleaning the Meter

Caring for the ACCU-CHEK Aviva meter is easy – just keep it free of dust. If you need to clean or disinfect it, follow these guidelines carefully to help you get the best performance possible:

Do

- Make sure the meter is off
- Gently wipe the meter's surface with a soft cloth slightly dampened (wring out any excess liquid) with one of these cleaning solutions:
 - Super Sani-Cloth
 - 70 % rubbing alcohol
 - Mild dishwashing liquid mixed with water
 - 10 % household bleach solution (1 part bleach plus 9 parts water) made the same day

Do Not

- Get any moisture in the code key slot or test strip slot
- Spray any cleaning solution directly onto the meter
- Put the meter under water or liquid
- Pour liquid into the meter

Maintenance and Testing

The meter needs little or no maintenance with normal use. It automatically tests its own systems every time you turn it on and lets you know if something is wrong. See Chapter 5, “Screen Messages and Troubleshooting.”

If you drop the meter or think it is not giving accurate results, make sure that your test strips and control solution haven't expired, then run a control test.

Test your display before each blood sugar test to make sure it is working properly. To test the display, turn off the meter, then press and **hold**  to see the complete display. All the segments should be clear and look exactly like the picture to the right. If one of the segments is missing or looks different from the picture, do not use the meter. Call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

Cleaning the Lancing Device

Clean the outside of the ACCU-CHEK Multiclix lancing device regularly with 70 % isopropyl (rubbing) alcohol. DO NOT place the entire device under water. DO NOT use bleach. At least once a week, disinfect the cap after cleaning by placing it in 70 % rubbing alcohol for 10 minutes. Allow the cap to air-dry after disinfecting.

Screen Messages and Troubleshooting


Never make treatment decisions based on an error message. If you have any concerns, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.


(blank screen)

The meter will not turn on or the display is blank.

- Battery is dead – Insert new battery.
- Display is damaged – Call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.
- Extreme temperatures – Move the meter to room temperature.

Battery power is low. Change the battery soon. See Chapter 5, “Changing the Battery.”

The meter is in the set-up mode, waiting for you to change or confirm settings. See Chapter 4, “Meter Setup.”

The meter is ready for you to insert a test strip.

Blood sugar may be higher than the measuring range of the system. See Chapter 3, “Understanding Test Results.”


The meter is ready for a drop of blood or control solution.


Blood sugar may be lower than the measuring range of the system. See Chapter 3, “Understanding Test Results.”

This test result was flagged. See Chapter 3, “Flagging Test Results” for more details.

The meter is not coded or the code key is not inserted. Turn off the meter and recode it. See Chapter 1, “Coding the Meter.”

Blood sugar is below the defined hypo (low blood sugar) level. See Chapter 3, “Understanding Test Results.”

The test strips expire at the end of the current month. At the end of the month, discard the code key and any remaining test strips. Insert a new code key from a new box of test strips and ensure the code key number

matches the code number on the test strip container. Make sure the time and date in the meter are correct.

Your blood sugar may be extremely low, or the test strip may be damaged or not properly inserted. If you see this error message **after** you applied blood to the test strip, see Chapter 3, "Unusual Test Results." If you see this error

message **before** you applied blood to the test strip, remove the test strip and reinsert it, or replace it if damaged. Verify that the code number printed on the container of test strips matches the code number printed on the code key currently inserted in the meter. If the message reappears, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

The code key is incorrect. Turn off the meter and insert a new code key. If this does not fix the problem, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

Your blood sugar may be extremely high or a meter or a test strip error has occurred. If this confirms how you feel, contact your doctor immediately. If it does not confirm how you feel, repeat

the test and refer to Chapter 3, “Unusual Test Results.” If this does not confirm how you feel, run a control test with your control solution and a new test strip. If the control result is within the acceptable range, review the proper testing procedure and repeat your blood sugar test with a new test strip. If the E-3 code still appears for your blood sugar test, your blood sugar result may be extremely high and above the system’s reading range. **Contact your doctor immediately.** If the control result is not within the acceptable range, see Chapter 2, “Understanding Control Test Results.”

Not enough blood or control solution was drawn into the test strip for measurement or was applied after the test has started. Discard the test strip and repeat the test.

The code key is from an expired lot of test strips. Ensure the code number matches the code number on the test strip container. Make sure the time and date in the meter are correct.

Blood or control solution was applied to the test strip before the flashing blood drop symbol appeared on the display. Discard the test strip and repeat the test.

The temperature is above or below the proper range for the meter. Refer to the test strip package insert for system operating conditions. Move to an area with the appropriate conditions, wait five minutes, and repeat the test. Do not artificially heat or cool the meter.

An electronic error occurred or, in rare cases, a used test strip was removed and reinserted. Turn the meter off and on or take the battery out for 20 seconds and reinsert it. Perform a blood sugar or control test. If the problem

persists, call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

The battery is almost out of power. Change the battery now. See Chapter 5, "Changing the Battery."

The time and date settings may be incorrect. Make sure the time and date are correct and adjust, if necessary. See Chapter 4, “Setting the Time and Date.”

Note: If you see any other error screen, please call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

Chapter 6: Technical Information

Product Limitations

Please read the literature packaged with the test strips to find the latest information on product specifications and limitations.

Specifications

Blood volume	See test strip insert for more information
Sample type	See test strip insert for more information
Measuring time	See test strip insert for more information
Measuring range	See test strip insert for more information
Test strip storage conditions	See test strip insert for more information
Meter storage conditions	Temperature: -13 °F to 158 °F
System operating conditions	See test strip insert for more information
Relative humidity operating range	See test strip insert for more information
Memory capacity	500 blood sugar results and 20 control results with time and date
Automatic power off	2 minutes
Power supply	One 3-volt lithium battery (coin cell 2032)
Display	LCD
Dimensions	94 x 53 x 22 mm (LWH)
Weight	Approximately 60 g (with battery)
Construction	Hand-held

Protection class	III
Meter type	The ACCU-CHEK Aviva meter is suitable for continuous operation
Lancing device depths	0.80, 0.95, 1.10, 1.25, 1.40, 1.55, 1.70, 1.85, 2.00, 2.15, 2.30 mm

Electromagnetic Compatibility

This meter meets the electromagnetic immunity requirements as per EN ISO 15197 Annex A. The chosen basis for electrostatic discharge immunity testing was basic standard IEC 61000-4-2. In addition, it meets the electromagnetic emissions requirements as per EN 61326. Its electromagnetic emission is thus low. Interference from other electrically driven equipment is not to be anticipated.

Performance Analysis

The performance data for the ACCU-CHEK Aviva system (ACCU-CHEK Aviva meter with ACCU-CHEK Aviva test strips) were obtained using capillary blood from diabetic patients (method comparison, accuracy), venous blood (repeatability) and control solution (reproducibility). The system is calibrated with venous blood containing various levels of glucose. The reference values are obtained using the hexokinase method. For method comparison, the results were compared with results obtained using the hexokinase method with deproteinization (automatic analyzer). The hexokinase method is traceable to an NIST standard.

The ACCU-CHEK Aviva system meets the EN ISO 15197 requirements.

Measuring Principle

Blood glucose concentrations may be measured in whole blood or plasma. Although you always apply whole blood to the test strip, the meter displays blood glucose results that relate to plasma. Please refer to the test strip package insert for information on how the system works, on the test principle, and on reference methods.

Product Safety Information

Strong electromagnetic fields may interfere with the proper operation of the meter. Do not use this meter close to sources of strong electromagnetic radiation.

To avoid electrostatic discharge, do not use the meter in a very dry environment, especially one in which synthetic materials are present.

Explanation of Symbols

CAUTION Please refer to safety-related notes in the manual accompanying this instrument.

3-volt battery coin cell 2032

The ACCU-CHEK Aviva meter has been listed by Underwriter's Laboratories, Inc.©, in accordance with UL 3101-1 and CAN/CSA C22.2 No. 1010-1

Store at

This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.

Travel Documentation

If you are traveling on a commercial airline, you may be required to provide documentation certifying that this meter meets environmental conditions and test procedures for Airborne Equipment (RTCA/DO-160E) section 21 Emission of Radio Frequency Energy. Please visit www.accu-chek.com or call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072 to obtain a copy of the document.

Disposing of the Meter, Test Strips, Lancets, and Batteries

Any product coming in contact with blood is considered contaminated (potentially infectious).* During normal testing, any blood glucose meter may come in contact with blood. Lancing devices may also be considered sharps. Disposal of sharps is regulated by law in many jurisdictions.

The European Union has a requirement for improving waste management practices for certain electronic equipment, but meters fall outside the scope of the European Directive 2002/96/EC.** This is not a requirement for the United States; however, Roche is committed to recycling and sustainability. Please consider the following points when disposing of your used testing materials:

- Comply with any laws or ordinances relating to the disposal of sharps and/or contaminated products. Contact your local health department or other appropriate authorities for proper handling and disposal of used meters, used test strips, used lancets, and used batteries.
- Consider recycling of the meters and batteries at an appropriate facility. Be aware the meter is potentially hazardous electronics scrap (e-scrap) and should be disposed of accordingly. The batteries are potentially hazardous also and should be disposed of accordingly.

- Decontaminate the meter before recycling or disposing. Wipe the outside of the meter with dilution of bleach solution (one part bleach to nine parts water).
- Users in professional environments (i.e., healthcare professionals) should follow the existing policies and procedures governing the proper handling and disposal of potentially infectious waste, e-scrap and batteries.

*29 CFR 1910.1030 – Bloodborne pathogens

**Directive 2002/96/EC – Directive on waste electrical and electronic equipment (WEEE)

Warranty

ACCU-CHEK Aviva Meter 30-day Money-back Guarantee

ACCU-CHEK Aviva Meter 30-day Money-back Guarantee for Qualifying Consumers

Roche Diagnostics offers qualifying consumers that purchase an ACCU-CHEK Aviva meter, a 30-day money back guarantee. If you are not fully satisfied with your ACCU-CHEK Aviva meter, call the ACCU-CHEK Customer Care Service Center toll-free at 1-800-858-8072 to determine whether you qualify to receive a full refund within 30 days of purchase. If you are covered under Medicare, Medicaid, other federal/state programs, or private insurance you are NOT eligible for this money-back offer. Consumers affected by this exclusion may instead request a different ACCU-CHEK meter/system. The refund will be limited to the amount paid by you net of any rebates. You must have a copy of the dated itemized purchase receipt and the original packaging to obtain this refund.

ACCU-CHEK Aviva Meter Limited 3-Year Warranty

Roche Diagnostics warrants to the original purchaser of the meter that your ACCU-CHEK Aviva meter will be free from defects in materials and workmanship for three years from the date of purchase. If, during this three-year period, the meter does not work properly because of a defect in materials or workmanship, Roche Diagnostics will replace it with a new ACCU-CHEK Aviva meter or equivalent product free of charge. The warranty on the replacement meter will expire on the date of the original warranty expiration or 90 days after the shipment of a replacement system, whichever period is longer. The purchaser's exclusive remedy with respect to the ACCU-CHEK Aviva meter shall be replacement.

This warranty does not apply to the performance of an ACCU-CHEK Aviva meter that has been damaged by accident or has been altered, misused, tampered with, or abused in any way. Roche Diagnostics will handle meters that show damage or abuse according to its Non-Warranty Service Policy described on the following page.

THE ABOVE WARRANTY IS EXCLUSIVE OF ALL OTHER WARRANTIES, AND ROCHE DIAGNOSTICS MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE TO THE PURCHASER OR ANY OTHER PERSON FOR ANY INCIDENTAL, CONSEQUENTIAL, INDIRECT, SPECIAL, OR PUNITIVE DAMAGES ARISING FROM OR IN ANY WAY CONNECTED WITH THE PURCHASE OR OPERATION OF THE METER OR ITS PARTS. NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IF ANY IS IMPLIED FROM THE SALE OF THE METER, SHALL EXTEND FOR A LONGER DURATION THAN THREE YEARS FROM THE DATE OF PURCHASE.

Some states do not allow limitations on how long an implied warranty will last or the exclusion of incidental or consequential damages, so the above limitation and exclusion may not apply to you. This warranty gives you specific legal rights, which vary from state to state.

Non-Warranty Service Policy

Roche Diagnostics Non-Warranty Service Policy applies to meters where the above warranty has not become effective, has become inapplicable, or has expired. Roche Diagnostics will replace, at its option, meters returned to it for a service charge (not to exceed \$35). Replacement will be with the same or similar product. Replacement meters will be warranted for a period of 90 days from shipment under a limited warranty providing for replacement of parts and labor at no charge.

Warranty and Service Instructions

All requests for return of ACCU-CHEK Aviva meters under the above warranty or service policy must be made to the ACCU-CHEK Customer Care Service Center. You will be mailed a return authorization label, which must be affixed to your carton for shipping the system to Roche Diagnostics. Cartons received without this label will be returned to you at your expense.

Customers experiencing difficulties should review the troubleshooting information in Chapter 5 of this booklet. Further inquiries should be directed to the ACCU-CHEK Customer Care Service Center.

Be sure to fill out and mail the Warranty Card that comes with the ACCU-CHEK Aviva system.

Limited License

CAUTION – A RESTRICTED LICENSE LIMITS USE OF THE ACCU-CHEK® AVIVA SYSTEM – READ CAREFULLY THE LIMITATIONS RECITED BELOW.

The ACCU-CHEK Aviva system (meter, including the code key, and test strips) and its use are protected by U.S. Patent Nos. 6,645,368; 5,366,609; 5,352,351; 5,122,244; 7,276,146; 7,276,147; 7,407,811; 7,452,457; 7,488,601; 7,494,816; and 7,569,126. A license to use the ACCU-CHEK Aviva system is only granted when the ACCU-CHEK Aviva meter is used with the ACCU-CHEK Aviva test strips and ACCU-CHEK Aviva code keys.

ACCU-CHEK Aviva test strips are provided with a specifically matched ACCU-CHEK Aviva code key. These test strips and code keys are specifically manufactured for operation with the ACCU-CHEK Aviva meter. Use of other test strips with an unmatched code key or even with the matched code key supplied by another manufacturer may prevent or impair the proper function of the ACCU-CHEK Aviva system.

Using the ACCU-CHEK Aviva system indicates your acceptance of the restricted license to use the ACCU-CHEK Aviva system only with ACCU-CHEK Aviva test strips and the ACCU-CHEK Aviva code key. Further, if you have purchased an ACCU-CHEK Aviva monitoring kit or an ACCU-CHEK Aviva meter that includes this restricted license, then this restricted license applies regardless of any additional offers found in ACCU-CHEK Aviva test strip packages. If you do not agree to the terms and conditions of the restricted license, you may return, at the place of purchase, the unused ACCU-CHEK Aviva system for a full refund. If you have any questions, please call the ACCU-CHEK Customer Care Service Center at 1-800-858-8072.

Except where prohibited by statute, all warranties covering the ACCU-CHEK Aviva system are voided by use of the ACCU-CHEK Aviva system with any test strips or code keys other than ACCU-CHEK Aviva test strips or code keys.

Patent Information

Aspects of the ACCU-CHEK Aviva system, including the meter, code key, and test strips, and their use are covered by one or more of the following United States patents: 5,352,351; 5,997,817; 5,438,271; 5,366,609; 6,645,368; 6,662,439; 5,122,244; 7,073,246; 7,276,146; 7,276,147; 7,338,639; 7,386,937; 7,407,811; 7,452,457; 7,488,601; 7,494,816; 7,556,723; 7,569,126.

CAUTION – A RESTRICTED LICENSE LIMITS USE OF THE ACCU-CHEK® MULTICLIX DEVICE (lancing device and lancets). READ CAREFULLY THE LIMITATIONS RECITED BELOW.

The ACCU-CHEK Multiclix device and its use are protected by U.S. Patent Nos. Re. 35,803; 6,419,661; 7,077,828; 7,223,276; 7,273,484; 7,322,998. Additional U.S. patents pending. A license to use the ACCU-CHEK Multiclix device is only granted when ACCU-CHEK Multiclix lancet drums are used with the ACCU-CHEK Multiclix device.

ACCU-CHEK Multiclix lancets are high precision components that are produced to the close tolerances required for satisfactory operation with the ACCU-CHEK Multiclix device. Use of other lancets with the ACCU-CHEK Multiclix device may prevent or impair proper function of the ACCU-CHEK Multiclix device.

Using the ACCU-CHEK Multiclix device indicates your acceptance of the restricted license to use the ACCU-CHEK Multiclix device only with ACCU-CHEK Multiclix lancet drums. Further, if you have purchased an

ACCU-CHEK Multiclix device that includes this restricted license, then this restricted license applies regardless of any additional offers found in ACCU-CHEK Multiclix device packages. If you do not agree to the terms and conditions of the restricted license, you may return, at the place of purchase, the unused ACCU-CHEK Multiclix device for a full refund. If you have any questions, please call the ACCU-CHEK Customer Care Services Center at 1-800-858-8072.

Except where prohibited by statute, all warranties covering the ACCU-CHEK Multiclix device are voided by use of the ACCU-CHEK Multiclix device with any lancet drums other than ACCU-CHEK Multiclix lancet drums.

Patent Information

The ACCU-CHEK Multiclix device and its use are covered by the following U.S. patents: Re. 35,803; 6,419,661; 7,077,828; 7,223,276; 7,273,484; 7,322,998.

Additional Supplies

The following supplies and accessories are available from your authorized Roche Diabetes Healthcare Center, pharmacies, or your medical/surgical supply dealer:

Test Strips

ACCU-CHEK Aviva 50 ct test strips

ACCU-CHEK Aviva 100 ct test strips

Control Solutions

ACCU-CHEK Aviva 2 level control solutions

Lancets

ACCU-CHEK Multiclix lancets 102 ct

Information for Healthcare Professionals

Healthcare Professionals: Follow the infection control procedures appropriate for your facility.

Sample Handling

Always wear gloves and follow your facility's infection control procedures when handling blood-contaminated items. Always adhere to the recognized procedures for handling objects that are potentially contaminated with human material. Follow the hygiene and safety policy of your laboratory or institution. Prepare the selected blood collection site per facility policy. A drop of fresh, whole blood is required to perform a blood glucose test. Fresh venous, capillary, arterial, or neonatal blood may be used. Caution should be taken to clear arterial lines before blood is drawn. Caution should be exercised in the interpretation of neonate blood glucose values below 50 mg/dL. Blood glucose determination with venous or arterial blood must be performed within 30 minutes of sample collection. Avoid air bubbles with the use of pipettes. For fresh venous blood specimens, refer to the test strip package insert for acceptable anticoagulants.

Recommending Alternative Site Testing to Patients

Decisions about whether to recommend alternative site testing (AST) should take into account the motivation and knowledge level of the patient and his or her ability to understand the considerations relative to diabetes and AST. If you are considering recommending AST for your patients, you need to

understand that there is a potential for a significant difference between fingertip and alternative site blood glucose test results. The difference in capillary bed concentration and blood perfusion throughout the body can lead to sample site-to-site differences in glucose results. These physiological effects vary between individuals and can vary within a single individual based upon his or her behavior and relative physical condition. Our studies involving AST of adults with diabetes show that most persons will find their glucose level changes more quickly in the fingers' blood than the alternative sites' blood.* This is especially important when glucose levels are falling or rising rapidly. If your patient is used to making treatment decisions based upon fingertip readings, he or she should consider the delay, or lag-time, affecting the reading obtained from an alternative site.

*Data on file

Index

- appetite, increased, 50
- battery, changing, 71
- battery, type, 72, 83, 85
- beeper settings, 59
- blood sugar, testing, 27
- blurred vision, 50
- button, on/off/set, 9
- code key, 9, 11
- comfort dial, 10, 30
- computer, connecting your meter to, 69
- control solution, 18
- control test results, understanding, 23
- control test, acceptable range, 23
- control test, results, 23
- control test, running, 19
- control test, unacceptable range, 23
- dehydration, 5
- display check, 74
- error messages, 78–81
- fatigue, 50
- frequent urination, 5, 50
- healthcare professionals, 95
- high blood sugar, 50
- hyperglycemia, 50
- hypoglycemia, 5, 50
- lancet drum, 10, 27
- limited license, 91–93
- low blood sugar, 50
- maintenance, meter, 74
- memory, meter, 53
- numbness, 50
- patent information, 91–93
- product limitations, 83
- product specifications, 83–84
- rapid heartbeat, 50
- screen messages, 75–81
- settings, meter, 56
- supplies, 94
- sweating, 50
- symbols, 85
- technical information, 83–87
- test results, understanding, 49
- test results, unusual, 49
- test strips, ACCU-CHEK Aviva, 10, 15
- thirst, increased, 50

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ACCU-CHEK, ACCU-CHEK AVIVA, ACCU-CHEK MULTICLIX, and SOFTCLIX are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.
©2010 Roche Diagnostics. All rights reserved.